Suprapubic Catheter for Single Use

Version: A/1

Before using the product, please read all information carefully

[Intended use]

The product is suitable for suprapubic cystostomy and urine drainage for clinical urinary diversion.

[Specification]

REF	Tip	Size	Balloon (ml)	Length (mm)	Material	Lumen
ED-270803	Central hole	8Fr	Integrated	3	310	Silicone 2-way
ED-271003	Central hole	10Fr	Integrated	3	310	Silicone 2-way
ED-271205	Central hole	12Fr	Integrated	5	400	Silicone 2-way
ED-271405	Central hole	14Fr	Integrated	5	400	Silicone 2-way
ED-271610	Central hole	16Fr	Integrated	10	400	Silicone 2-way
ED-271810	Central hole	18Fr	Integrated	10	400	Silicone 2-way
ED-272010	Central hole	20Fr	Integrated	10	400	Silicone 2-way
ED-272210	Central hole	22Fr	Integrated	10	400	Silicone 2-way
ED-272410	Central hole	24Fr	Integrated	10	400	Silicone 2-way

[Direction for use]

- Preparation: Remove catheter from sterile pouch using aseptic technique. Inflate the balloon with a syringe without need and check the leakage of the balloon and the valve. Deflate the balloon completely.
- The bladder should be completely filled up.
- Carry out an incision of the skin.
- Remove the protective tube of the puncture cannula.
- Insert the catheter up to the first (distal) marker, however only as far until the tip of the catheter is positioned within the tip of the puncture cannula.
- Introduce the cannula with inserted catheter through the skin incision vertical to the abdominal wall until the bladder is punctured and urine drains off.
- Advance the catheter through the puncture cannula and hold the puncture cannula firmly until the second marker has reached the cannula to ensure that the balloon is completely inside the bladder.
- Connect the catheter to a suitable urinary drainage bag.
- Pull out the puncture cannula above skin level.
- Inflate the balloon with a syringe without needle with sterile aqua dest. with 10% glycerine solution
- To remove the catheter, place the syringe firmly on the valve and deflate the balloon completely. Pull the catheter out carefully.

[Contraindication]

Haemorrhagic diathesis, venous congestion, varicosis, open or infectious skin injuries around the puncture site, empty bladder, ascites, peritoneal carcinomas, lower abdominal tumours, bladder cancer, pregnancy.

[Precaution]

- The cutting edge and tip of the puncture cannula must not damage the catheter.
- Safely dispose of the puncture cannula. Connect catheter to the urine drainage bag.
- This product has been sterilized by ethylene oxide gas, and is intended for single use only.
- If packing is damaged, do not use.

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UT_IFU-SUL-EN_Rev2.0 (11-23)

- The size and balloon capacity is marked on the outer unit package and the funnel of the catheter.
- When placing suprapubic catheter, the catheter may not be withdrawn through the puncture cannula..
- The scalpel-like ends of the puncture cannula might destroy the catheter tip.
- This product should be operated by medical staff.
- Indwelling time recommend no more than 30 days.

[Symbol]**[Warning]**

When inflate the balloon, shall not exceed the nominal capacity printed on the catheter (ml).

[Storage]

Store in a cool, dark and dry place, the temperature should not be higher than 40 C with good ventilation, but without corrosive gas.

[Date of language] The product is only used in EU countries

[Date of manufacture] See inner packing label

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CH REP Obelis Swiss GmbH
Ruessenstrasse 12, 6340 Baar/ZG, Switzerland

UK REP MedPath Limited
27 Old Gloucester Street, London, United Kingdom, WC1N 3AX

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Urotech GmbH
Produttore: HAIYAN KANGYUAN MEDICAL INSTRUMENT CO.,LTD.
EC REP Lotus NL B.V.